


MedSource
I.V. Infusionw

ClearSafe **Comfort**[™]
Safety I.V. Catheter



The ClearSafe Comfort[™] Safety I.V. Catheter was designed to allow care givers easy and safe I.V. access with cost effective solutions. The ClearSafe Comfort[™] I.V. Catheter is manufactured with PUR for maximum patient comfort. Our specially designed beveled needle ensures smooth, painless insertion. The clicking lock mechanism ensures safety from cross contamination. MedSource has simplified needle stick safety and compliance.

Item #	Description	UOM
MS-84214	ClearSafe Comfort I.V. Catheter, 14g x 1 1/4"	200/Case
MS-84216	ClearSafe Comfort I.V. Catheter, 16g x 1 1/4"	200/Case
MS-84218	ClearSafe Comfort I.V. Catheter, 18g x 1 1/4"	200/Case
MS-84220	ClearSafe Comfort I.V. Catheter, 20g x 1 1/4"	200/Case
MS-842202	ClearSafe Comfort I.V. Catheter, 20g x 1"	200/Case
MS-84222	ClearSafe Comfort I.V. Catheter, 22g x 1"	200/Case
MS-84224	ClearSafe Comfort I.V. Catheter, 24g x 3/4"	200/Case

Features:

- Fully encapsulating needle shield and clicking lock safety mechanism
- Ergonomically designed thumb grip
- Gently tapered catheter tip
- Large-volume chamber

Benefits:

- Designed to help healthcare workers as well as hospital waste handlers from exposure to blood borne pathogens via needlestick injuries.
- Allows for a variety of user techniques to aid insertion.
- Provides smooth transition from needle to catheter.
- Allows clear visual conformation through flashback and easy handling.

Specification: ClearSafe Comfort™ Safety I.V. Catheter

- Manufactured with PUR (Polyurethane)
- Siliconised Ground Stainless Steel Needle

- Smooth surface for easy catheter insertion and advancement.
- Softer, more comfortable for the patient and longer in-dwelling performance.
- Sharp needle technology provides both ease of penetration and improved patient comfort during cannulation.