



MEDIPAK UVLI PROTECTION

- UV Light Inhibiting Film In Green & Amber
- Designed to Block 90% or Greater UV Light
- Engineered for High Clarity
- Regular UVLI Covers for Transport & Short-Term Storage
- Slit-Top UVLI Covers for Hanging IV Bags
- UVLI Tubing Covers for Line Protection
- UVLI Zip Bags for Resealable Transport
- Tamper-Evident UVLI Bags to Support Chain of Custody/Diversion Reduction
- Multiple Sizes Available

